KNEE PROSTHESIS SYSTEM INSTRUCTIONS FOR USE (IFU) (EN)

CE₂₈₀₃

EN

IMPORTANT NOTE:

The HIPKNEE Knee system (Proknee-Orthoknee) is a reconstructive knee system intended for patients suffering from pain and/or deformity of the knee joint.

DESCRIPTION OF THE SYSTEM:

HIPKNEE manufactures this Knee System in a variety of sizes to accommodate differing surgeon and patient requirements.

The Knee System consists of :

- Femoral components (Cemented and Uncemented) manufactured from cast Cobalt Chromium alloy confirming to ISO 5832-4 - ASTM F 75
- Tibial components (Cemented and Uncemented) manufactured from cast Cobalt Chromium alloy confirming to ISO 5832-4 - ASTM F 75
- Tibial Component Stem manufactured from Titanium Alloy confirming to ISO 5832-3
- Tibial inserts manufactured from UHMWPE / Cross Linked UHMWPE confirming ISO 5834-2
- Patella components manufactured from UHMWPE confirming ISO 5834-2

Description of the component material is provided on the labels. Each Knee Component is designed as part of a system and does not allow the substitution of components from other systems or manufacturers. All implantable devices are designed sterile and for Single-Use only. The Hydroxyapatite or Titanium Plasma / Porous coatings applied to implant surfaces are intended for uncemented arthroplasty.

Intended Use: The HIPKNEE Knee Prosthesis System is intended for use in knee arthroplasty to reduce pain, correct deformity and improve knee function and mobility in skeletally mature patients with damaged knee joint surfaces.

Intended User: The primary user of this device is an orthopaedic surgeon trained and experienced in knee arthroplasty. The user shall have adequate professional qualification, training and experience in total knee arthroplasty, revision knee arthroplasty, cemented/cementless fixation techniques and related surgical instrumentation.

Indications:

Knee Systems are indicated for use in knee arthroplasty for reduction or relief of pain and/or improved knee function in skeletally mature patients with the following conditions:

- 1) Non-inflammatory degenerative joint disease including osteoarthritis, traumatic arthritis, or avascular necrosis;
- 2) inflammatory degenerative joint disease including rheumatoid arthritis;
- 3) correction of functional deformity;

Indications For Cemented And Cementless Fixation:

Cemented Prostheses: Cemented femoral components and cemented tibial components are intended for use where immediate fixation is achieved through the use of bone cement (PMMA). Cemented fixation may be preferred in patients with reduced bone quality, osteoporosis, poor bone stock, elderly patients, or where the surgeon determines that cemented fixation provides the most appropriate clinical outcome.

Cementless Prostheses: Cementless femoral components and cementless tibial components are intended for biological fixation through bone ongrowth and/or ingrowth. Cementless components are recommended in patients with adequate bone quality and sufficient bone stock to achieve primary stability and long-term biological fixation.

Hydroxyapatite (HA), Titanium Plasma or Porous Coated Components:

Components incorporating Hydroxyapatite (HA), Titanium Plasma or porous coatings are designed for cementless fixation and shall not be implanted using bone cement unless specifically indicated on the product label.

The surgeon is responsible for selecting the appropriate fixation method, cemented or cementless, based on patient anatomy, bone quality, age, activity level, clinical condition and intraoperative findings.

POSSIBLE ADVERSE EFFECTS:

The general principles of good patient selection and sound surgical judgment apply to the total knee procedure. Preoperative planning and meticulous surgical technique are essential to achieve optimum results. Considerations of anatomic loading, soft-tissue condition, and component placement are critical to minimize a variety of postoperative complications.

1. Osteolysis (progressive bone resorption). Osteolysis can be asymptomatic, and therefore, routine periodic radiographic examination is vital to prevent any serious future complication
2. Particulate generation leading to increased wear rates necessitating early revision. Soft tissue imbalance leading to excessive wear;
3. Allergic reactions to materials; metal sensitivity that may lead to histological reactions;
4. Delayed wound healing; deep wound infection (early or late) which may necessitate removal of the prosthesis. On rare occasions, arthrodesis of the involved joint or amputation of the limb may be required;
5. A sudden drop in blood pressure intra-operatively due to the use of bone cement;
6. Damage to blood vessels or hematoma;
7. Temporary or permanent nerve damage, peripheral neuropathies and subclinical nerve damage as possible result of surgical trauma resulting in pain or numbness of the affected limb;
8. Cardiovascular disorders including venous thrombosis, pulmonary embolism, or myocardial infarction;
9. Dislocation, migration and/or subluxation of prosthetic components from improper positioning, trauma, loss of fixation and/or muscle and fibrous tissue laxity;
10. Periarticular calcification or ossification, with or without impediment to joint mobility;
11. Varus-valgus deformity;
12. Traumatic arthrosis of the knee from intraoperative positioning of the extremity;
13. Inadequate range of motion due to improper selection or positioning of components, periarticular calcification, flexion contracture;
14. Femoral, tibial or patellar bone or component fracture intraoperatively or postoperatively; fracture by trauma or excessive loading, particularly in the presence of poor bone stock;
15. Undesirable shortening or lengthening of the limb;

16. Aggravated problems of the affected limb or contralateral extremity by leg length discrepancy, excess femoral medialization, or muscle deficiency;
17. Pain.

Contraindications

1. Tibia, femur or patella fractures.
2. Blood Circulatory constriction such as arteriosclerosis
3. Known allergic reaction to CoCr
4. Infection
5. Children or patients with skeletal immaturity
6. Tumors, cysts in the distal femur, proximal tibia
7. Cementless components are contraindicated in patients with insufficient bone stock or bone quality where adequate primary stability cannot be achieved. Cemented components are contraindicated in patients with known hypersensitivity to bone cement constituents.

WARNINGS AND PRECAUTIONS:

The patient should be warned of surgical risks, and made aware of possible adverse effects. The patient should be warned that the devices do not replace normal healthy bone, and that the implant can break or become damaged as a result of strenuous activity or trauma, and has a finite expected service life and may need to be replaced in the future. Higher rates of wear may be initiated by particles of cement, metal, or other debris which can cause abrasion of the articulating surfaces. Higher rates of wear may shorten the useful life of the prosthesis, and lead to early revision surgery to replace the worn prosthetic components. Loosening, bending, cracking, or fracture of implant components; fracture of the implant can occur as a result of trauma, strenuous activity.

The potential long-term biological effects of metal wear debris and metal ion production are not known. Questions regarding carcinogenicity have been raised in literature; no studies have conclusive evidence that metal wear debris or metal ions are carcinogenic.

Ligament stability, based on maintained isometry throughout knee of motion. Designed to maintain collateral stability throughout the range of motion. Improper selection, placement, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may lead to subsequent reduction in the service life of the prosthetic components. Malalignment of the components or inaccurate implantation can lead to excessive wear and/or failure of the implant or procedure. Inadequate preclosure cleaning (removal of surgical debris) can lead to excessive wear. Use clean gloves when handling implants. Laboratory testing indicates that implants subjected to body fluids, surgical debris or fatty tissue have lower adhesion strength to cement than implants handled with clean gloves. Improper preoperative or intraoperative implant handling or damage (scratches, dents, etc.) can lead to crevice corrosion, fretting, fatigue fracture and/or excessive wear. Do not modify implants. The surgeon is to be thoroughly familiar with the surgical technique, implants and instruments prior to performing surgery.

MR Information

The HIPKNEE Knee Prosthesis System is MR Conditional. Before MRI examination, the patient should inform the physician and MRI personnel about the implant. MRI examination should be performed only when the conditions for

an MR Conditional implant are met. The implant may affect MR image quality around the implant area.

Preoperative:

1. Patient selection should consider the following factors which could lead to increased risk of failure and can be critical to the eventual success of the procedure: the patient's weight, activity level, occupation and lifestyle.
2. The patient should be warned of surgical risks, and made aware of possible adverse effects. The patient should be warned that the prosthesis does not replace normal healthy bone, that the prosthesis can break or become damaged as a result of certain activity or trauma, has a finite expected service life, and may need to be replaced at some time in the future. The patient should also be advised of other risks that the surgeon believes should be disclosed. The patient should be advised that any noise or unusual sensation should be reported to the surgeon as it may indicate implant malfunction.
3. Use care in handling and storing of implant components. Cutting, bending, or scratching the surfaces of components can significantly reduce the strength, fatigue resistance and/or wear characteristics of the implant system. These in turn may induce internal stresses that are not obvious to the eye and may lead to fracture of the component.
4. Surgical information is available upon request. The surgeon should be familiar with the technique.
5. An adequate inventory of implant sizes should be available at the time of surgery.
6. Intraoperative fracture or breaking of instruments can occur. Instruments which have experienced extensive use or excessive force are susceptible to fracture. Instruments should be examined for wear and damage prior to surgery.

Intraoperative:

1. The correct selection of the implant is extremely important. The appropriate type and size should be selected for patients with consideration of anatomical and biomechanical factors such as patient age and activity levels, weight, bone and muscle conditions, etc. Generally, the largest cross-section component which will allow adequate bone support to be maintained is preferred. Failure to use the optimum size component may result in loosening, bending, cracking, or fracture of the component and/or bone.
2. Modular components must be assembled securely to prevent disassociation. Avoid repeated assembly and disassembly of the modular components which may compromise a critical locking action of the components. Surgical debris must be cleaned from components before assembly. Debris may inhibit the proper fit and locking of modular components leading to early failure of the procedure.
3. Prior to closure, the surgical site should be thoroughly cleaned of bone chips, extraneous cement, etc.; foreign particles at the metal and/or plastic interface may cause excessive wear and/or friction.
4. An implant should never be reused. While it may appear undamaged, imperfection may exist which would reduce the service life of the implant.

Extramedullary alignment metot:

Rotational alignment can be achieved by aligning the ankle clamp attachment shaft with second toe. You can choose 5 degree or 0 degree posterior angle on resection guide.

Postoperative:

1. The patient must be advised of the limitations of the reconstruction and the need for protection of the prosthesis from full weight bearing until adequate fixation and healing have occurred. Excessive activity and trauma affecting the joint replacement have been implicated with failure of the reconstruction by loosening, fracture and/or wear of the prosthetic components. Loosening of the components can result in increased production of wear particles, as well as damage to the bone, making successful revision surgery more difficult.
2. Periodic, long-term follow-up is recommended to monitor the position and state of the prosthetic components, as well as the condition of the adjoining bone.
3. Periodic postoperative x-rays are recommended for close comparison with early post-op conditions to detect long term evidence of changes in position, loosening, bending, or cracking of components

PACKAGING AND LABELING:

Components should only be accepted if received by the hospital or surgeon with the factory packaging and labeling intact. Products cannot be reused as indicated on the label.

STERILIZATION / RESTERILIZATION:

All implants are supplied sterile and have been packaged in protective trays. The method of sterilization is noted on the package label. All radiation sterilized components have been exposed to a minimum of 25 kiloGrays of gamma radiation. Inspect packages for punctures or other damage prior to surgery. If the sterile barrier has been broken, return the component to HIPKNEE. If metal products are out of use after opening, they cannot be sterilized again. HIPKNEE accepts no responsibility for the use of implants resterilized after contact with body tissues or fluids. Polyethylene products cannot be resterilized.

MEANING OF SYMBOLS:

	Consult the user manual.
	Single sterile barrier system with protective packaging outside
	Lot numbers
	Product code number (Catalog number)
	Manufacturer name and address
	Sterilized by Gamma Irradiation
	Use on time (Expiration date)
	Not reusable

	Production date
	European Union symbol of conformity
	Medical device
	Indicates a carrier that contains unique device identifier information
	Upper limit of temperature
	Do not resterilize
	Do not use if package is damaged and consult instructions for use
	Caution
	Humidity limitation
	European Article Number
	MR Conditional

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